Diamond TexturedNitrile Gloves (Blue)



PRODUCT INFORMATION					
MATERIAL	Nitrile				
COLOR	Blue				
ТҮРЕ	Ambidextrous, non-sterile, single-use				
INTERIOR	Powder-free				
EXTERIOR	Fully diamond textured				
SIZES	8 (M), 9(L), 10(XL), 11(2XL)				
COUNTRY OF ORIGIN	Thailand				
STORAGE	Store in original packaging in a cool, dry and well ventilated area, away from dust, direct sunlight, moisture, x-ray and excessive heat above 100°F (37°C)				



PHYSICAL PROPERTIES					
AQL	1.5				
GLOVE WEIGHT	14.1g (size 8)				
GLOVE THICKNESS	9mil				
GLOVE LENGTH	12"				
	BEFORE AGING	AFTER AGING			
TENSILE STRENGTH (MPA)	min. 22	min. 22			
ULTIMATE ELONGATION	min. 500% min. 400%				





QUALITY STANDARDS			
FDA STATUS	(21 CFR 177) compliant for food handling		
AUDIT STANDARDS	Manufactured in an ISO 9001:2015 and an ISO 13485:2016 facility		
TEST STANDARDS	EN ISO 374-1:2016/Type B EN ISO 374-5:2016 Resistance to Bacteria, Fungi & Virus		

PACKAGING & ORDERING INFORMATION					
CODE	SIZE	PURCHASE UNIT	CASE DIMENSIONS (LxWxH)	CASE WEIGHT	CUBIC FEET
102130	8		15.5 X 10.5 X 12"	19.4lbs	0.49ft ³
102140	9	1 case of 500 Gloves			
102150	10	(50/bag x 10)			
102160	11				

RESISTANCE OF GLOVES TO PERMEATION BY CHEMICALS									
CHEMICAL			EN ISO 374-1:2016 PERFORMANCE LEVEL		EN 374-4:2019 MEAN DEGRADATION / %				
Sodium Hydroxide 40% (K)			6		7.3				
Sulphuric Acid 96% (L)			1		100.0				
Ammonium Hydroxide 25% (O)			2		6.1				
Formaldehyde 37% (T)			6		0.9				
Hydrogen Peroxide 30% (P)			2		10.5				
EN ISO 374-1:2016 - permeation levels are based on breakthrough times as follows:									
Performance Level:	1	2	3	4	5	6			
Minimum breakthrough time (Min):	>10	>30	>60	>120	>240	>480			

Safety gloves to protect against chemicals are classified according to their permeation time (time taken for the chemical to penetrate the glove) and number of chemicals tested:

- Type A at least 30min each for at least 6 test chemicals
- Type B at least 30min each for at least 3 test chemicals
- Type C at least 10min each for at least 1 test chemicals

EN 374-4:2019 - Degradation results indicate the change in puncture resistance of the gloves after exposure to the challenge chemical

EN ISO 374-5:2016 - Resistance to Bacteria and Fungi = Pass, Resistance to Virus = Pass

MANDATORY STATEMENTS EN ISO 374-1:2016

"This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals."

"The chemical resistance has been assessed under laboratory conditions from samples taken from the palm only (except in cases where the glove is equal to or over 400mm - where the cuff is tested also) and relates only to the chemical tested. It can be different if the chemical is used in a mixture."

"It is recommended to check that the gloves are suitable for the intended used because the conditions at the workplace may differ from the type depending on temperature, abrasion and degradation."

"When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves."

"The penetration resistance has been assessed under laboratory conditions and relates to the tested specimen."



Contact us today to receive samples or for more information on this product.

Eagle Protect PBC

Uncle Mats and Supply 954-751-9800 www.UncleSupply.com info@UncleSupply.com





SAFEGUARD WHAT MATTERS TH